

Measuring Adsorption Processes on Platinum Electrodes using the Solartron ModuLab Potentiostat.

John Harper
Solartron Analytical, Farnborough, UK.

Application Guide: ANML01

Introduction

The determination of the active surface area of Platinum is of great interest to electrochemists. For example fuel cell scientists wish to calculate the utilization of catalyst to improve performance and reduce costs. This is traditionally achieved using cyclic voltammetry as the diagnostic technique. Furthermore, the principles presented in this guide can be applied to the analysis of any cyclic voltammogram. This demonstration guide explains how such calculations are performed using the ModuLab software package.

Experimental

i) Chemicals and Electrodes

- WE: 2 mm Platinum Disk (Goodfellows >99.9 %) (working area = 0.0314 cm²), CE: Platinum mesh, RE: Hg / Hg₂SO₄
- 0.5 M H₂SO₄

ii) Hardware Required

- 1 M sample s⁻¹ ModuLab Potentiostat

iii) Software Settings

- 1st Vertex Potential = 1.2 V vs. Ref, 2nd Vertex Potential = -0.35 V vs. Ref
- Scan rate = 100 mV s⁻¹
- Measurement Mode: E change (amount of change = 5 mV)

Results and Discussion

A sample voltammogram is presented in Figure 1 with the oxygen and hydrogen regions identified. The shaded area represents the region associated with the desorption of hydrogen including double layer charging. The total charge is calculated using the following expression;

$$Q = \frac{1}{\nu} \int_{E_1}^{E_2} I dE \quad (1)$$

where Q = charge (C), ν = scan rate (V s⁻¹), E₁ = lower potential limit, E₂ = upper potential limit, I = current (amps)

To perform this calculation in ModuLab, follow these steps;

- Select the region of the curve to integrate. The software requires two markers to define this region; marker A defines E₁ and B defines E₂ (see Figure 1).
- (See Figure 2) change graph 1 to Q vs. t in the axes control panel and click 'Refresh'. Note that the cursors are at the same point with respect to time as those in Figure 1.
- To calculate the total charge, subtract the value of the charge at point A from that at point B.

(Note: These values will appear on the screen when you hover over the green text boxes).

This is the total charge associated with this region of the I-V curve. One must now subtract the double layer charge.




Figure 1: Cyclic Voltammogram of a Pt disk electrode in H₂SO₄. The hydrogen desorption region is highlighted in the shaded region.

Figure 3 shows the expanded I-V curve of the region of interest. The dashed regions (Cdl (a) and Cdl (b)) represent the double layer charging whilst the region H is charge associated with the desorption of hydrogen. To calculate Cdl (a) and Cdl (b);

- i) first define the baseline. Using the line fit tool in the fitting options, place two cursors on a linear section of the curve and click 'run'. The results of the regression analysis will be presented in the dialogue box together with the line fit in the graph.
- ii) Using the equation of a straight line, $y = mx + c$, calculate the value of I_2 at $x = E_2$. Similarly, calculate I_1 at $x = E_1$
- iii) For the area Cdl (b),

$$Q = \frac{1}{2\nu} \int_{E_1}^{E_2} I.dE \text{ or } [(I_2 - I_1) * (E_2 - E_1)] / 2 * \nu$$

- iv) Similarly, the area Cdl (a) ,

$$Q = \frac{1}{\nu} \int_{E_1}^{E_2} I.dE \text{ or } [(I_1) * (E_2 - E_1)] \nu$$

- v) Subtract both of these values from the total charge to yield the charge associated with the desorption step.




Figure 2: Plot of Q vs t (derived from Figure 1 in the software)

Worked Example

Using the data from the results in this note, the following calculation is intended to guide the user through the methodology presented in the previous sections.

i) Using Figure 2 to calculate the total charge: $A = 9.32 \times 10^{-5} \text{ C}$, $B = 1.18 \times 10^{-4} \text{ C}$.
Therefore total charge $(B - A) = 24.8 \mu\text{C}$

ii) Using figure 3 to calculate double layer charging:

Results of regression analysis, $(y = mx + c)$

$$m = -1.58 \times 10^{-6}, c = 1.134 \times 10^{-6}$$

Therefore substituting for x, we can calculate the following

$$I_2 = -1.58 \times 10^{-6} \times E_{21} (0.34 \text{ V}) + 1.134 \times 10^{-6} = 1.7 \mu\text{A}$$

$$I_1 = 1.1 \mu\text{A} \text{ (at } E_1 = 0\text{V)}$$

iii) $\text{Area Cdl (b)} = [(I_2 - I_1) * |E_2 - E_1|] / 2$
 $= [(1.7 \mu\text{A}) - (1.1 \mu\text{A}) * (0.34)] / 2 * 0.1$
 $= 1.02 \mu\text{C}$

iv) $\text{Area Cdl (a)} = [I_1 * |E_2 - E_1|]$
 $= (1.1 \mu\text{A}) * (0.34) / 0.1$
 $= 3.7 \mu\text{C}$

Charge due to desorption = $24.8 \mu\text{C} - 4.7 \mu\text{C} = 20.1 \mu\text{C}$ for a 2 mm diameter electrode

*** or $639 \mu\text{C cm}^{-2}$. ***

According to the literature, the charge associated with the monolayer adsorption of hydrogen = $210 \mu\text{C cm}^{-2}$. The difference in results reflect the difference between the geometric and true surface area of the electrode under investigation.




Figure 3: Expanded I-V curve of hydrogen desorption region. (derived from Figure 1). The total charge, Q, comprises the charge associated with desorption (region H) and that of double layer charging (Cdl (a) and (b)).

Conclusions

The ModuLab software package allows the rapid and accurate determination of adsorption/desorption processes on electrode surfaces. This will be of great benefit to fundamental and applied electrochemical studies.



Solartron Analytical's Quality System is approved to BS EN ISO 9001:1994



FM01709

...part of **AMETEK**® Advanced Measurement Technology

UNIT B1 ARMSTRONG MALL
SOUTHWOOD BUSINESS PARK
FARNBOROUGH, GU14 0NR
UNITED KINGDOM
Phone: +44 (0) 1252 556 800
Fax: +44 (0) 1252 556 899

801 SOUTH ILLINOIS AVENUE
OAK RIDGE
TN 37831-2011
USA
Phone: +1 865 425 1360
Fax: +1 865 425 1334

Visit our website for a complete list of our global offices and authorized agents

solartron.info@ametek.com

www.solartronanalytical.com